Wide Range of External Interfacing

Connects seamlessly with administration PCs or external systems. Also provides easy connectivity with barcode readers, card readers or other convenient devices for loading patient data for examinations.



TX-20P Specifications

| Туре | Non-contact |
|------------------------|---|
| Non-contact tonometer | |
| Measurement range | 0 to 60 mmHg (automatic switching between 30/60 mmHg) Units: mmHg/hPa |
| Measuring increment | 1 mmHg |
| Eye fixation target | Internal fixation lamp / External fixation lamps |
| Non-contact pachymeter | |
| Measurement range | 400 to 1300 µm |
| Measuring increment | 10 µm |
| Measurement mode | Full Auto / Auto / Manual |
| | Full Auto / Auto: 3D Auto alignment |
| Operation distance | 11 mm (from front of nozzle) |

| Data input/output | Input: USB |
|-------------------|---|
| | Output: LAN/RS-232C (D-Sub 9-pin male) |
| Display | 5.7-Inch color LCD monitor |
| Printer | Built-in thermal line printer with auto cutter |
| Power supply | AC100 to 240 V 50/60 Hz 0.8 to 0.4 A |
| Operating range | Front/back: 40 mm, Left/right: 90 mm, Up/down: 30 mm, Chin rest up/down: 60 mm |
| Dimensions | W260 mm × D490 mm × H500 mm |
| Weight | 15 kg |
| | |

| | Main Components |
|---|------------------------------|
| [| TX-20P main unit |
| | Nozzle cap |
| | Printing paper |
| | Chin rest paper (100 sheets) |
| Γ | Dust cover |

Options

Chin rest paper (500 sheets)

Canon

 Manufacturer
 CANON INC.
 9-1, Imaikami-cho, Nakahara-ku, Kawasaki, Kanagawa 211-8501, Japan
 Telephone: (81)-3-3758-2111



Make sure you read the manual before using the instrument.
Use the correct power source and voltage shown.
Make sure the instrument is properly connected to ground. Failure to do so may result in an electric shock if there is a malfunction or ground leakage.

* Specifications and appearance are subject to change without notice.

Canon



Full Auto continuous measurements of both eyes One-touch operation for measuring intraocular pressure and corneal thickness

Compact body, with functions for simultaneously measuring intraocular pressure and corneal thickness. Designed for ease-of-use with tiltable color LCD monitor and other convenient features. TX-20P, the full auto non-contact tonometer that is hygienic for both patient and examiner.



Intraocular Pressure and Corneal Thickness Measurements with a Single Unit

In addition to intraocular pressure measurement, the TX-20P also includes a function for measuring central corneal thickness with the use of a split light system. The corneal thickness measurement line is shown based on the corneal thickness measurement results. This helps prevent underestimation of intraocular pressure in patients with thinning cornea such as after having had LASIK.

Designed for Smooth Measurements

Full Auto mode coupled with simple and intuitive controls makes for smooth measurements, even for those using the unit for the first time.

Advanced intelligent 3D alignment

In Full Auto mode, a single tap on the measurement button at the start of the examination quickly and easily measures the intraocular pressure and corneal thickness of both the left and right eyes, and automatically prints out the results. In addition to Full Auto, Auto and Manual measurement modes can be selected.

Compact, Lightweight and Energy Efficient

At just 15 kg, the weight of the TX-20P main unit has been reduced significantly. Its small maximum area of lateral movement translates into excellent space efficiency, and it is also energy efficient.

Diverse Functions Designed for Usability

Color LCD monitor with tilt function

Equipped with a 5.7-inch color LCD monitor that can be swiveled by up to around 40° for the ease of view for the examiner. This allows for stress-free measurements to be taken, even while the examiner is standing.



Auto Cut printout

A new Auto Cut function has been added to the TX-20P's built-in printer. The printout slot at the front also boosts convenience.



Corneal thickness measurement function

Equipped with a central corneal thickness measurement function (pachymeter) based on a slit light system. This allows both the intraocular pressure and central corneal thickness to be measured quickly and easily at the same time. Measurements of intraocular pressure based on corneal thickness mean it is an effective method for glaucoma screening and other ailments.



Equipped with external eye fixation lamps

In addition to the internal fixation lamp, the TX-20P also features external fixation lamps which can be switched easily depending on the condition of the patient's eye.



Motorized joystick

The TX-20P is equipped with a motorized joystick with measurement button, making operation much easier.



Intraocular pressure measurement function

The TX-20P measures intraocular pressure and central corneal thickness at the same time, and uses the obtained results to automatically calculate and display the compensated intraocular pressure. This helps prevent underestimation of intraocular pressure in patient eyes with thinner cornea due to having had LASIK, or overestimation of intraocular pressure in patient eyes with thicker corneas.

Corneal thickness compensated intraocular pressure 9.3 Standard corneal 538 540 Most recent corneal 0 Menu 10 Menu 10

nickness